

For Research Use Only

SSBP4 Polyclonal antibody

Catalog Number: 15595-1-AP



Basic Information

Catalog Number: 15595-1-AP	GenBank Accession Number: BC000274	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 µg/ml by Nanodrop and 347 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 170463	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: single stranded DNA binding protein 4	
Isotype: IgG	Calculated MW: 39 kDa	
Immunogen Catalog Number: AG7952	Observed MW: 39-42 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : U-937 cells,
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



U-937 cells were subjected to SDS PAGE followed by western blot with 15595-1-AP (SSBP4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.