

For Research Use Only

KCNH6 Polyclonal antibody

Catalog Number: 15640-1-AP



Basic Information

Catalog Number: 15640-1-AP	GenBank Accession Number: BC006334	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 120 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 81033	
Source: Rabbit	Full Name: potassium voltage-gated channel, subfamily H (eag-related), member 6	
Isotype: IgG	Calculated MW: 994aa, 110 kDa; 502aa, 56 kDa	
Immunogen Catalog Number: AG8107		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data