

For Research Use Only

# COL5A1 Polyclonal antibody

Catalog Number: 15657-1-AP

1 Publications



## Basic Information

Catalog Number:

15657-1-AP

Size:

150ul, Concentration: 160 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8009

GenBank Accession Number:

BC008760

GeneID (NCBI):

1289

UNIPROT ID:

P20908

Full Name:

collagen, type V, alpha 1

Calculated MW:

184 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

## Notable Publications

| Author          | Pubmed ID | Journal             | Application |
|-----------------|-----------|---------------------|-------------|
| Chengjuan Zhang | 39755111  | Int Immunopharmacol | WB          |

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

