

For Research Use Only

GRB7 Polyclonal antibody

Catalog Number: 15738-1-AP



Basic Information

Catalog Number: 15738-1-AP	GenBank Accession Number: BC006535	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 167 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2886	
Source: Rabbit	Full Name: growth factor receptor-bound protein 7	
Isotype: IgG	Calculated MW: 532 aa, 60 kDa	
Immunogen Catalog Number: AG8537		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Grb7 belongs to a family of adapter proteins. This family, including Grb10, Grb7 and Grb 14, contains a proline-rich, a RA, a PH, a BPS and SH2 domain. Grb7 associate with a number of receptors including EGFR and FAK. Functions of Grb7 include mobility because it associates receptors such FAK, and direct evidence from overexpression of this proteins. This antibody recognize the C-terminal of GRB7.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data