

For Research Use Only

LOC81691 Polyclonal antibody

Catalog Number: 15762-1-AP



Basic Information

Catalog Number:

15762-1-AP

Size:

150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8393

GenBank Accession Number:

BC007646

GeneID (NCBI):

81691

UNIPROT ID:

Q96IC2

Full Name:

exonuclease NEF-sp

Calculated MW:

774 aa, 87 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human kidney tissue were subjected to SDS PAGE followed by western blot with 15762-1-AP (LOC81691 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.