

For Research Use Only

DDX24 Polyclonal antibody, PBS Only

Catalog Number:15769-1-PBS

Featured Product



Basic Information

Catalog Number:

15769-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8407

GenBank Accession Number:

BC008847

GeneID (NCBI):

57062

UNIPROT ID:

Q9GZR7

Full Name:

DEAD (Asp-Glu-Ala-Asp) box

polypeptide 24

Calculated MW:

96 kDa

Observed MW:

120 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

DDX24 is a family member of Asp-Glu-Ala-Asp (DEAD) box containing RNA helicases, and DEAD-box RNA helicases are characterized by a conserved DEAD motif. DDX24 was associated with cancer development, viral infection, and vascular malformation.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




HT-29 cells were subjected to SDS PAGE followed by western blot with 15769-1-AP (DDX24 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15769-1-PBS in a different storage buffer formulation.




Various lysates were subjected to SDS PAGE followed by western blot with 15769-1-AP (DDX24 antibody) at dilution of 1:7000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15769-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using DDX24 antibody (15769-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red). This data was developed using the same antibody clone with 15769-1-PBS in a different storage buffer formulation.



WB result of DDX24 antibody (15769-1-AP; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-DDX24 transfected A549 cells. This data was developed using the same antibody clone with 15769-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 15769-1-AP (DDX24 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 15769-1-PBS in a different storage buffer formulation.