

For Research Use Only

GLRX Polyclonal antibody, PBS Only

Catalog Number: 15804-1-PBS

Featured Product



Basic Information

Catalog Number:

15804-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8535

GenBank Accession Number:

BC005304

GeneID (NCBI):

2745

UNIPROT ID:

P35754

Full Name:

glutaredoxin (thioltransferase)

Calculated MW:

106 aa, 12 kDa

Observed MW:

12 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15804-1-AP (GLRX antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15804-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed HepG2 cells using GLRX antibody (15804-1-AP) at dilution of 1:200 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 15804-1-PBS in a different storage buffer formulation.

Various lysates were subjected to SDS PAGE followed by western blot with 15804-1-AP (GLRX antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15804-1-PBS in a different storage buffer formulation.