

For Research Use Only

SGTB Polyclonal antibody

Catalog Number: 15895-1-AP



Basic Information

Catalog Number:

15895-1-AP

Size:

150ul, Concentration: 160 ug/ml by Nanodrop and 133 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8709

GenBank Accession Number:

BC012044

GeneID (NCBI):

54557

UNIPROT ID:

Q96EQ0

Full Name:

small glutamine-rich tetratricopeptide repeat (TPR)-containing, beta

Calculated MW:

304 aa, 33 kDa

Observed MW:

33-36 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15895-1-AP (SGTB antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.