

For Research Use Only

IMPA2 Polyclonal antibody

Catalog Number: 15925-1-AP



Basic Information

Catalog Number:

15925-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8755

GenBank Accession Number:

BC017176

GeneID (NCBI):

3613

UNIPROT ID:

O14732

Full Name:

inositol(myo)-1(or 4)-
monophosphatase 2

Calculated MW:

288 aa, 31 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

