

For Research Use Only

QRSL1 Polyclonal antibody, PBS Only

Catalog Number: 15926-1-PBS



Basic Information

Catalog Number:

15926-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8756

GenBank Accession Number:

BC014389

GeneID (NCBI):

55278

UNIPROT ID:

Q9HOR6

Full Name:

glutamyl-tRNA synthase (glutamine-hydrolyzing)-like 1

Calculated MW:

528 aa, 57 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 15926-1-AP (QRSL1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15926-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 15926-1-AP (QRSL1 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 15926-1-PBS in a different storage buffer formulation.