

For Research Use Only

MRPS18C Polyclonal antibody

Catalog Number: 15928-1-AP



Basic Information

Catalog Number:

15928-1-AP

Size:

150ul , Concentration: 300 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8771

GenBank Accession Number:

BC005186

GeneID (NCBI):

51023

UNIPROT ID:

Q9Y3D5

Full Name:

mitochondrial ribosomal protein

S18C

Calculated MW:

142 aa, 16 kDa

Observed MW:

16 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP : HepG2 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-MRPS18C (IP:15928-1-AP, 4ug;
Detection:15928-1-AP 1:300) with HepG2 cells
lysate 1925 ug.