

For Research Use Only

# UBE2Q1 Polyclonal antibody

Catalog Number: 16072-1-AP



## Basic Information

**Catalog Number:**

16072-1-AP

**Size:**

150ul , Concentration: 650 ug/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG9028

**GenBank Accession Number:**

BC009286

**GeneID (NCBI):**

55585

**UNIPROT ID:**

Q7Z7E8

**Full Name:**

ubiquitin-conjugating enzyme E2Q family member 1

**Calculated MW:**

422 aa, 46 kDa

**Observed MW:**

46 kDa

**Purification Method:**

Antigen affinity purification

**Recommended Dilutions:**

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

## Applications

**Tested Applications:**

IP, ELISA

**Species Specificity:**

human

**Positive Controls:**

IP : HEK-293 cells,

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



IP result of anti-UBE2Q1 (IP:16072-1-AP, 4ug;  
Detection:16072-1-AP 1:600) with HEK-293 cells  
lysate 1440 ug.