

For Research Use Only

STRA13 Polyclonal antibody

Catalog Number: 16116-1-AP



Basic Information

Catalog Number:

16116-1-AP

Size:

150ul , Concentration: 20 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9044

GenBank Accession Number:

BC009571

GeneID (NCBI):

201254

UNIPROT ID:

A8MT69

Full Name:

stimulated by retinoic acid 13 homolog (mouse)

Calculated MW:

63 aa, 7 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

STRA13, also named as CENPX, FAAP10 and MHF2, is a centromere protein with MW 10 kDa. It is recruited to forks stalled by DNA interstrand cross-links, and required for cellular resistance to such lesions. STRA13 may be involved as a repressor in a number of decision events occurring during differentiation of various cell lineages.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

