

For Research Use Only

PHYHD1 Polyclonal antibody, PBS Only

Catalog Number: 16125-1-PBS



Basic Information

Catalog Number:

16125-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9108

GenBank Accession Number:

BC009373

GeneID (NCBI):

254295

UNIPROT ID:

Q5SRE7

Full Name:

phytanoyl-CoA dioxygenase domain containing 1

Calculated MW:

297 aa, 33 kDa

Observed MW:

30-32 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 16125-1-AP (PHYHD1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16125-1-PBS in a different storage buffer formulation.