

For Research Use Only

ZMIZ2-Specific Polyclonal antibody

Catalog Number: 16158-1-AP



Basic Information

Catalog Number: 16158-1-AP	GenBank Accession Number: NM_031449	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 83637	
Source: Rabbit	Full Name: zinc finger, MIZ-type containing 2	
Isotype: IgG	Calculated MW: 97 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

ZMIZ2 and ZMIZ1 members of a PIAS-like family of proteins that interact with nuclear hormone receptors. This antibody is specific for ZMIZ2.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data