

For Research Use Only

# CXCR4B-specific Polyclonal antibody



Catalog Number: 16161-1-AP

## Basic Information

<b>Catalog Number:</b> 16161-1-AP	<b>GenBank Accession Number:</b> NM_001008540	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 7852	
<b>Source:</b> Rabbit	<b>Full Name:</b> chemokine (C-X-C motif) receptor 4	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 40 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

CXCR4B, also named as CXCR4-LO, CD184, LESTR, LCR1, FB22, NPYRL, HM89 and Fusin. The antibody is specific to CXCR4B isoform.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data