

For Research Use Only

MAGEC1-Specific Polyclonal antibody



Catalog Number: 16176-1-AP

Basic Information

Catalog Number:

16176-1-AP

Size:

150ul , Concentration: 473 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005462

GeneID (NCBI):

9947

UNIPROT ID:

O60732

Full Name:

melanoma antigen family C, 1

Calculated MW:

124 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

This gene is a member of the melanoma antigen gene (MAGE) family. The proteins of this family are tumor-specific antigens. Expression of this protein was not seen in normal tissues, with the exception of testis. This protein contains a large number of unique short repetitive sequences in front of the MAGE-homologous sequence, and therefore is about 800 aa longer than the other MAGE proteins. (From OMIM)

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

