

For Research Use Only

# GMEB1-Specific Polyclonal antibody

Catalog Number: 16189-1-AP



## Basic Information

**Catalog Number:**

16189-1-AP

**Size:**

150ul , Concentration: 33 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_006582

**GeneID (NCBI):**

10691

**UNIPROT ID:**

Q9Y692

**Full Name:**

glucocorticoid modulatory element binding protein 1

**Calculated MW:**

63 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human

## Background Information

GMEB1 is a trans-acting factor that binds to glucocorticoid modulatory elements (GME), acting in estrogen receptor-regulated pathway GREB1 is expressed in proliferating prostate tissue & prostate cancer and is thought to function in hormone-dependent cancer growth. It associates with GMEB2 protein, and the complex is essential for parvovirus DNA replication. The protein acts complicatedly in SDS-PAGE and may have multiple electrophoresis forms.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

