

For Research Use Only

GMEB1-Specific Polyclonal antibody

Catalog Number: 16189-1-AP



Basic Information

Catalog Number: 16189-1-AP	GenBank Accession Number: NM_006582	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 33 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 10691	
Source: Rabbit	Full Name: glucocorticoid modulatory element binding protein 1	
Isotype: IgG	Calculated MW: 63 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

GMEB1 is a trans-acting factor that binds to glucocorticoid modulatory elements (GME), acting in estrogen receptor-regulated pathway GREB1 is expressed in proliferating prostate tissue & prostate cancer and is thought to function in hormone-dependent cancer growth. It associates with GMEB2 protein, and the complex is essential for parvovirus DNA replication. The protein acts complicatedly in SDS-PAGE and may have multiple electrophoresis forms.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data