

For Research Use Only

Collagen Type XVIII Polyclonal antibody



Catalog Number: 16191-1-AP

1 Publications

Basic Information

Catalog Number:

16191-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_030582

GeneID (NCBI):

80781

UNIPROT ID:

P39060

Full Name:

collagen, type XVIII, alpha 1

Calculated MW:

178 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF, IHC

Species Specificity:

human

Cited Species:

human, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Hao Zhou	38515609	Bioact Mater	IHC, IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

