

For Research Use Only

FOXS1 Polyclonal antibody

Catalog Number: 16234-1-AP



Basic Information

Catalog Number: 16234-1-AP	GenBank Accession Number: BC013408	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 453 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2307	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: forkhead box S1	
Isotype: IgG	Calculated MW: 330 aa, 35 kDa	
Immunogen Catalog Number: AG8957	Observed MW: 35 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : mouse liver tissue, rat liver tissue
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 16234-1-AP (FOXS1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.