

For Research Use Only

CEP27 Polyclonal antibody

Catalog Number: 16298-1-AP



Basic Information

Catalog Number:

16298-1-AP

Size:

150ul , Concentration: 187 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8898

GenBank Accession Number:

BC010903

GeneID (NCBI):

55142

UNIPROT ID:

Q9NVX0

Full Name:

centrosomal protein 27kDa

Calculated MW:

235 aa, 27 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

CEP27, also named as C15orf25 and HAUS2, is a centrosomal protein with MW 27 kDa. It contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

