## For Research Use Only

## ELSPBP1 Polyclonal antibody

Catalog Number: 16316-1-AP



**Basic Information** 

Catalog Number:

GenBank Accession Number:

**Purification Method:** 

Antigen affinity purification

16316-1-AP

BC015598 GeneID (NCBI):

150ul , Concentration: 67  $\mu g/ml$  by Bradford method using BSA as the

64100 **UNIPROT ID:** Q96BH3

standard; Source:

Full Name: epididymal sperm binding protein 1

Rabbit Isotype:

Calculated MW: 223 aa, 26 kDa

Immunogen Catalog Number:

**Tested Applications:** 

**ELISA** 

Species Specificity:

human

Storage

**Applications** 

Storage:

Store at -20°C. Stable for one year after shipment.

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

Selected Validation Data