

For Research Use Only

RAB30 Polyclonal antibody

Catalog Number: 16328-1-AP



Basic Information

Catalog Number: 16328-1-AP	GenBank Accession Number: BC014213	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 27314	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: RAB30, member RAS oncogene family	
Isotype: IgG	Calculated MW: 203 aa, 23 kDa	
Immunogen Catalog Number: AG9444	Observed MW: 20-23 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, PC-12 cells
Species Specificity: human, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 16328-1-AP (RAB30 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.