

For Research Use Only

HYI Polyclonal antibody

Catalog Number: 16390-1-AP



Basic Information

Catalog Number:

16390-1-AP

Size:

150ul , Concentration: 600 ug/ml by Nanodrop and 327 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9457

GenBank Accession Number:

BC019041

GeneID (NCBI):

81888

UNIPROT ID:

Q5T013

Full Name:

hydroxypyruvate isomerase homolog (E. coli)

Calculated MW:

217aa, 24 kDa; 174aa, 19 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse liver tissue, rat liver tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16390-1-AP (HYI antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.