

For Research Use Only

HYI Polyclonal antibody, PBS Only

Catalog Number:16390-1-PBS



Basic Information

Catalog Number:

16390-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9457

GenBank Accession Number:

BC019041

GeneID (NCBI):

81888

UNIPROT ID:

Q5T013

Full Name:

hydroxypyruvate isomerase homolog (E. coli)

Calculated MW:

217aa,24 kDa; 174aa,19 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16390-1-AP (HY1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16390-1-PBS in a different storage buffer formulation.