

For Research Use Only

# PI4KAP2 Polyclonal antibody

Catalog Number: 16422-1-AP **1 Publications**



## Basic Information

**Catalog Number:**

16422-1-AP

**Size:**

150ul , Concentration: 93 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG9619

**GenBank Accession Number:**

BC020225

**GeneID (NCBI):**

375133

**Full Name:**

phosphatidylinositol 4-kinase, catalytic, alpha pseudogene 2

**Calculated MW:**

539aa,61 kDa; 577aa,66 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Cited Applications:**

WB

**Species Specificity:**

human, mouse, rat

**Cited Species:**

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Safiyah Ziyad	29386109	Cell Rep	WB

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data