

For Research Use Only

RBM38 Polyclonal antibody, PBS Only

Catalog Number: 16435-1-PBS



Basic Information

Catalog Number:

16435-1-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9748

GenBank Accession Number:

BC018711

GeneID (NCBI):

55544

UNIPROT ID:

Q9H0Z9

Full Name:

RNA binding motif protein 38

Calculated MW:

216 aa, 23 kDa

Observed MW:

25-31 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16435-1-AP (RBM38 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16435-1-PBS in a different storage buffer formulation.