

For Research Use Only

GPLD1 Polyclonal antibody

Catalog Number: 16444-1-AP **1 Publications**



Basic Information

| | | |
|---|---|--|
| Catalog Number: 16444-1-AP | GenBank Accession Number: BC020748 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 80 µg/ml by Bradford method using BSA as the standard; | GeneID (NCBI): 2822 | |
| Source: Rabbit | Full Name: glycosylphosphatidylinositol specific phospholipase D1 | |
| Isotype: IgG | Calculated MW: 176aa,20 kDa; 840aa,92 kDa | |
| Immunogen Catalog Number: AG9782 | | |

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|------------|-----------|---------|-------------|
| Indu Sinha | 27027327 | OMICS | WB |

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data