

For Research Use Only

Caveolin-1 Polyclonal antibody

Catalog Number: 16446-1-AP **3 Publications**



Basic Information

Catalog Number: 16446-1-AP	GenBank Accession Number: BC006432	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 857	
Source: Rabbit	UNIPROT ID: Q03135	
Isotype: IgG	Full Name: caveolin 1, caveolae protein, 22kDa	
Immunogen Catalog Number: AG10924	Calculated MW: 22 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Xiansheng Ye	32167773	J Agric Food Chem	WB
Jorge L Jimenez-Macias	40009687	Sci Adv	WB
Xiaohui Sun	37716034	Phytomedicine	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data