

For Research Use Only

# SNAPC2 Polyclonal antibody

Catalog Number: 16506-1-AP



## Basic Information

<b>Catalog Number:</b> 16506-1-AP	<b>GenBank Accession Number:</b> BC011868	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 167 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 6618	
<b>Source:</b> Rabbit	<b>Full Name:</b> small nuclear RNA activating complex, polypeptide 2, 45kDa	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 334 aa, 36 kDa	
<b>Immunogen Catalog Number:</b> AG9689		

## Applications

**Tested Applications:**  
IF, ELISA

**Species Specificity:**  
human

## Background Information

SNAPC2, also named as SNAP45, is a subunit of SNAPc complex. It recruits TBP and BRF2.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data