

For Research Use Only

HSP10 Polyclonal antibody

Catalog Number: 16512-1-AP



Basic Information

Catalog Number: 16512-1-AP	GenBank Accession Number: BC023518	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 450 µg/ml by Nanodrop;	GeneID (NCBI): 3336	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: heat shock 10kDa protein 1 (chaperonin 10)	
Isotype: IgG	Calculated MW: 102 aa, 11 kDa	
Immunogen Catalog Number: AG9732	Observed MW: 11 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16512-1-AP (HSP10 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.