

For Research Use Only

HSP10 Polyclonal antibody

Catalog Number: 16512-1-AP **1 Publications**



Basic Information

Catalog Number:

16512-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9732

GenBank Accession Number:

BC023518

GeneID (NCBI):

3336

UNIPROT ID:

P61604

Full Name:

heat shock 10kDa protein 1 (chaperonin 10)

Calculated MW:

102 aa, 11 kDa

Observed MW:

11 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

chicken

Positive Controls:

WB : HeLa cells, HepG2

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16512-1-AP (HSP10 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.