

For Research Use Only

IMPACT Polyclonal antibody

Catalog Number: 16528-1-AP



Basic Information

Catalog Number:

16528-1-AP

Size:

150ul , Concentration: 300 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9776

GenBank Accession Number:

BC034016

GeneID (NCBI):

55364

UNIPROT ID:

Q9P2X3

Full Name:

Impact homolog (mouse)

Calculated MW:

320 aa, 36 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HEK-293 cells, TT cells, HL-60 cells, LNCaP cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16528-1-AP (IMPACT antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.