

For Research Use Only

MRPS9 Polyclonal antibody, PBS Only

Catalog Number: 16533-1-PBS

Featured Product



Basic Information

Catalog Number:

16533-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9809

GenBank Accession Number:

BC057240

GeneID (NCBI):

64965

UNIPROT ID:

P82933

Full Name:

mitochondrial ribosomal protein S9

Calculated MW:

396 aa, 46 kDa

Observed MW:

38-45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16533-1-AP (MRPS9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16533-1-PBS in a different storage buffer formulation.



WB result of MRPS9 antibody (16533-1-AP; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-MRPS9 transfected HEK-293 cells. This data was developed using the same antibody clone with 16533-1-PBS in a different storage buffer formulation.