

For Research Use Only

MRP4 Polyclonal antibody

Catalog Number: 16568-1-AP



Basic Information

Catalog Number:
16568-1-AP

Size:
150ul, Concentration: 153 µg/ml by
Bradford method using BSA as the
standard;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG9850

GenBank Accession Number:
BC041560

GeneID (NCBI):
10257

Full Name:
ATP-binding cassette, sub-family C
(CFTR/MRP), member 4

Calculated MW:
859 aa, 97 kDa

Purification Method:
Antigen affinity purification

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data