

For Research Use Only

NEURL3 Polyclonal antibody, PBS Only

Catalog Number: 16648-1-PBS



Basic Information

Catalog Number:

16648-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10035

GenBank Accession Number:

BC012317

GeneID (NCBI):

93082

UNIPROT ID:

Q96EH8

Full Name:

neuralized homolog 3 (Drosophila) pseudogene

Calculated MW:

167 aa, 18 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-NEURL3 (IP:16648-1-AP, 3ug;
Detection:16648-1-AP 1:300) with mouse brain
tissue lysate 5000ug. This data was developed
using the same antibody clone with 16648-1-PBS in
a different storage buffer formulation.