For Research Use Only

LSM12 Polyclonal antibody

Catalog Number: 16670-1-AP



Purification Method:

Antigen affinity purification

Basic Information

Catalog Number:

16670-1-AP

Size:

standard;

Source:

Isotype:

BC105288

GenBank Accession Number:

GeneID (NCBI):

150ul , Concentration: 67 $\mu g/ml$ by Bradford method using BSA as the

124801 **UNIPROT ID:**

Q3MHD2 Full Name:

Rabbit

LSM12 homolog (S. cerevisiae) Calculated MW:

131 aa, 14 kDa

Immunogen Catalog Number:

AG10150

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

This product is exclusively available under Proteintech Group brand and is not available to purchase from any

other manufacturer.

Selected Validation Data