

For Research Use Only

TCL1B Polyclonal antibody

Catalog Number: 16793-1-AP



Basic Information

Catalog Number:

16793-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10264

GenBank Accession Number:

BC041616

GeneID (NCBI):

9623

UNIPROT ID:

O95988

Full Name:

T-cell leukemia/lymphoma 1B

Calculated MW:

128 aa, 15 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

TCL1b is mainly restricted to the lymphoid cell lineage and to the early stages of embryogenesis where it co-activates PKB, possibly to promote a growth advantage during development through PKB-stimulated cell survival.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data