

For Research Use Only

TCL1B Polyclonal antibody

Catalog Number: 16793-1-AP



Basic Information

Catalog Number: 16793-1-AP	GenBank Accession Number: BC041616	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 9623	
Source: Rabbit	Full Name: T-cell leukemia/lymphoma 1B	
Isotype: IgG	Calculated MW: 128 aa, 15 kDa	
Immunogen Catalog Number: AG10264		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

TCL1b is mainly restricted to the lymphoid cell lineage and to the early stages of embryogenesis where it co-activates PKB, possibly to promote a growth advantage during development through PKB-stimulated cell survival.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data