

For Research Use Only

# M6PR Polyclonal antibody

Catalog Number: 16795-1-AP **2 Publications**



## Basic Information

<b>Catalog Number:</b> 16795-1-AP	<b>GenBank Accession Number:</b> BC024206	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 67 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 4074	
<b>Source:</b> Rabbit	<b>Full Name:</b> mannose-6-phosphate receptor (cation dependent)	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 277 aa, 31 kDa	
<b>Immunogen Catalog Number:</b> AG10267		

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
IF, WB

**Species Specificity:**  
human, mouse, rat

**Cited Species:**  
mouse, rat

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Anup Parchure	36173346	J Cell Biol	IF
Jian Wei	28658628	Cell Rep	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data