

For Research Use Only

RP9 Polyclonal antibody, PBS Only

Catalog Number: 16832-1-PBS



Basic Information

Catalog Number:

16832-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10441

GenBank Accession Number:

BC025928

GeneID (NCBI):

6100

UNIPROT ID:

Q8TA86

Full Name:

retinitis pigmentosa 9 (autosomal dominant)

Calculated MW:

221 aa, 26 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 16832-1-AP (RP9 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16832-1-PBS in a different storage buffer formulation.



WB result of RP9 antibody (16832-1-AP; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-RP9 transfected HepG2 cells. This data was developed using the same antibody clone with 16832-1-PBS in a different storage buffer formulation.