

For Research Use Only

DEFB129 Polyclonal antibody

Catalog Number: 16847-1-AP



Basic Information

Catalog Number:

16847-1-AP

Size:

150ul , Concentration: 247 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9943

GenBank Accession Number:

BC070359

GeneID (NCBI):

140881

UNIPROT ID:

Q9H1M3

Full Name:

defensin, beta 129

Calculated MW:

183 aa, 20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

DEFB129, also named as C20orf87 and DEFB29, has antibacterial activity.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

