

For Research Use Only

# RAG1AP1 Polyclonal antibody

Catalog Number: 16884-1-AP



## Basic Information

**Catalog Number:**

16884-1-AP

**Size:**

150ul , Concentration: 73 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG10350

**GenBank Accession Number:**

BC009621

**GeneID (NCBI):**

55974

**UNIPROT ID:**

Q9BRV3

**Full Name:**

recombination activating gene 1 activating protein 1

**Calculated MW:**

167aa, 19 kDa; 221aa, 25 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human, mouse, rat

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

