

For Research Use Only

# RAG1AP1 Polyclonal antibody

Catalog Number: 16884-1-AP



## Basic Information

**Catalog Number:**  
16884-1-AP

**Size:**  
150ul , Concentration: 73 µg/ml by  
Bradford method using BSA as the  
standard;

**Source:**  
Rabbit

**Isotype:**  
IgG

**Immunogen Catalog Number:**  
AG10350

**GenBank Accession Number:**  
BC009621

**GeneID (NCBI):**  
55974

**Full Name:**  
recombination activating gene 1  
activating protein 1

**Calculated MW:**  
167aa,19 kDa; 221aa,25 kDa

**Purification Method:**  
Antigen affinity purification

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data