

For Research Use Only

PKN3 Polyclonal antibody, PBS Only

Catalog Number: 16904-1-PBS



Basic Information

Catalog Number:

16904-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10518

GenBank Accession Number:

BC062558

GeneID (NCBI):

29941

UNIPROT ID:

Q6P5Z2

Full Name:

protein kinase N3

Calculated MW:

889 aa, 99 kDa

Observed MW:

99 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MDA-MB-231 cells were subjected to SDS PAGE followed by western blot with 16904-1-AP (PKN3 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16904-1-PBS in a different storage buffer formulation.