

For Research Use Only

ACTRT2 Polyclonal antibody

Catalog Number: 16992-1-AP

3 Publications



Basic Information

Catalog Number:

16992-1-AP

Size:

150ul , Concentration: 150 ug/ml by Nanodrop and 133 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10597

GenBank Accession Number:

BC029499

GeneID (NCBI):

140625

UNIPROT ID:

Q8TDY3

Full Name:

actin-related protein T2

Calculated MW:

376 aa, 42 kDa

Observed MW:

40-42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, CoIP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Xiao-Zhen Zhang	35616329	Development	WB
Xiao-Zhen Zhang	35793634	Cell Rep	CoIP
Qi Zhang	38414365	Hum Reprod	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16992-1-AP (ACTRT2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.