

For Research Use Only

RPL18 Polyclonal antibody, PBS Only

Catalog Number:17029-1-PBS



Basic Information

Catalog Number:

17029-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10474

GenBank Accession Number:

BC009708

GeneID (NCBI):

6141

UNIPROT ID:

Q07020

Full Name:

ribosomal protein L18

Calculated MW:

188 aa, 22 kDa

Observed MW:

22-26 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17029-1-AP (RPL18 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17029-1-PBS in a different storage buffer formulation.