

PRSS23 Polyclonal antibody

Catalog Number: 17031-1-AP

Basic Information

Catalog Number:

17031-1-AP

Size:

150ul , Concentration: 260 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10491

GenBank Accession Number:

BC001278

GeneID (NCBI):

11098

UNIPROT ID:

O95084

Full Name:

protease, serine, 23

Calculated MW:

383 aa, 43 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

PRSS23(Serine protease 23) is also named as ZSIG13(Putative secreted protein Zsig13) and belongs to the peptidase S1 family. Enhanced PRSS23 expression is observed in various types of cancers, including breast, prostate, Papillary thyroid, and pancreatic cancers, and the expression of the PRSS23 has been linked with tumor progression in human(PMID:22291950).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

