

For Research Use Only

# HBEGF Polyclonal antibody

Catalog Number: 17051-1-AP

Featured Product

1 Publications



## Basic Information

### Catalog Number:

17051-1-AP

### Size:

150ul , Concentration: 127 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG10716

### GenBank Accession Number:

BC033097

### GeneID (NCBI):

1839

### UNIPROT ID:

Q99075

### Full Name:

heparin-binding EGF-like growth factor

### Calculated MW:

208 aa, 23 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Species Specificity:

human, mouse, rat

### Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
X Wang	25321193	Br J Cancer	

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

